



Reviews

The Orsiro Ultrathin, Bioresorbable-Polymer Sirolimus-Eluting Stent: A Review of Current Evidence

Brian J. Forrestal, Brian C. Case, Charan Yerasi, Hector M. Garcia-Garcia, Ron Waksman *

Section of Interventional Cardiology, MedStar Washington Hospital Center, Washington, DC, United States of America

ARTICLE INFO

Article history:

Received 2 October 2019

Received in revised form 19 December 2019

Accepted 30 December 2019

Keywords:

Coronary artery disease

Biodegradable-polymer stent

Drug-eluting stent

Acute coronary syndrome

ABSTRACT

Advances in stent design and the development of bioresorbable polymers have allowed the development of novel stent technologies such as the Orsiro bioresorbable-polymer sirolimus eluting stent (BP-SES). Over several non-inferiority trials, the BP-SES has demonstrated itself to be a safe and effective therapy for obstructive coronary artery disease. This article reviews the current evidence of the efficacy of the BP-SES and examines its performance in high-risk populations, such as patients presenting with ST-segment myocardial infarction, chronic total occlusions, diabetes, and small vessel disease.

© 2020 Elsevier Inc. All rights reserved.

1. Introduction

In recent years, a series of advances in stent technology has led to steady improvements in rates of restenosis and stent thrombosis and has reduced the need for repeat revascularization. The development of ultrathin struts, improvements in polymer biocompatibility, and the development of bioresorbable polymers have all aimed to decrease early stent thrombosis and reduce local inflammatory and hypersensitivity reactions, leading to reduced rates of in-stent restenosis, need for revascularization, and lower incidence of other adverse clinical events [1]. Biodegradable-polymer stents were designed to allow the steady and controlled eluting of antiproliferative drugs in the early phases of treatment. This mitigates neointimal hyperplasia [2] while at the same time, through the polymer's degradation, removing the nidus of inflammation at the vessel interface that is thought to be responsible for the “catch up” phenomenon of very late stent thrombosis observed with early generation drug-eluting stents (DES) [3].

The Orsiro BP-SES is available in diameters of 2.25, 2.5, 2.75, 3.0, 3.5, and 4.0 mm, and is available in lengths ranging from 9 to 40 mm. Stents

with diameters of 2.25 to 3.0 mm have a strut thickness of 60 µm, while stents with diameters of 3.5 to 4.0 mm have a strut thickness of 80 µm [4]. Drawing together recent technological advances, the Orsiro coronary stent (Biotronik, Buelach, Switzerland) consists of two distinct layers packaged around its metallic struts [2].

The innermost layer of the system consists of a cobalt-chromium alloy PRO-Kinetic energy™ stent arranged in a double helix pattern, which is designed to improve flexibility and deliverability. The middle proBIO™ layer consists of a silicon carbide coating that seals the metal alloy surface and reduces the allergenic interaction of metal ions with the vessel wall and blood pool [5]. Finally, the outer BIOLute™ layer is composed of a bioabsorbable poly-L-lactic acid (PLLA) polymer containing the antiproliferative sirolimus agent. The active BIOLute™ coating is distributed in an asymmetric fashion on the struts, with a thickness of 7.5 µm on the abluminal side and a thinner 3.5-µm layer on the luminal portion of the stent. Fig. 1 demonstrates the stent strut design. The PLLA polymer undergoes a hydrolytic reaction upon contact with the blood pool, causing the polymer to break down into CO₂ and H₂O via the Krebs cycle. The sirolimus load is 1.4 µg/mm² [6] with the elution optimized for 12 to 14 weeks after implantation. Fig. 2 summarizes some of the key design features in comparison with other contemporary stent designs.

In vitro studies have shown that within 30 days, 50% of the drug is eluted, and within 3 months, 80% is eluted. The PLLA degrades over 2 years, leaving the thin-strut cobalt-chromium stent with minimal endothelial injury [5]. This hybrid approach of combining ultrathin struts, a passive intermediate layer, and bioabsorbable polymer is aimed at creating less flow disturbance, limiting the long-term inflammatory

Abbreviations: BP-SES, biodegradable-polymer sirolimus-eluting stent; CTO, chronic total occlusion; CI, confidence interval; DCB, drug-coated balloon; DES, drug-eluting stent; DP-DES, durable-polymer drug-eluting stent; LLL, late lumen loss; PLLA, poly-L-lactic acid polymer; PP-EES, permanent-polymer everolimus-eluting stent; STEMI, ST-elevation myocardial infarction; TLF, target lesion failure; TLR, target lesion revascularization; TVF, target vessel failure; TVR, target vessel revascularization.

* Corresponding author at: MedStar Washington Hospital Center, 110 Irving St., NW, Suite 4B-1, Washington, DC 20010, United States of America.

E-mail address: ron.waksman@medstar.net (R. Waksman).

response at the vessel interface and encouraging rapid endothelialization with early vascular healing [1,2,5]. In larger stent diameters, where the comparative advantage of strut thickness is less (80 μ m in stents with diameters >3.5 mm), the additional design features of the inert silicon carbide layer and biodegradable polymer may still confer benefits in vascular healing and long-term inflammation in comparison with other contemporary stents of similar strut thickness.

2. Discussion

Several investigator- and industry-initiated randomized controlled trials have examined the efficacy and safety of the Orsiro BP-SES in recent years. Cumulatively, these trials have enrolled close to 30,000 patients around the world and have examined de novo lesions, acute coronary syndromes, and high-risk subgroups such as small-vessel disease, diabetes, and chronic total occlusions. Table 1 describes the key randomized trials and registries, which are discussed below, followed by an additional focused discussion on the Orsiro BP-SES performance in high-risk subgroups.

2.1. The BIOFLOW trials

The introduction of drug-eluting stents (DES) brought about major advances in combating early restenosis but also raised concerns about delayed intimal healing and prolonged inflammation [7] at the site of implant, leading to late events such as stent thrombosis and in-stent restenosis [8]. The development of biodegradable-polymer stents was intended to reduce the local inflammatory response and allow improved endothelialization, decreasing the risk of late events such as stent thrombosis and in-stent restenosis [9]. Starting in 2013, the BIOFLOW trials were conducted to examine the safety, efficacy, and noninferiority of the Orsiro BP-SES against several leading second-generation DES. Results of these studies were used to support regulatory approval for the use of the Orsiro stent system in Europe, Japan, China, and the United States.

BIOFLOW I [6] was a first-in-man study designed to evaluate the safety and effectiveness of the newly developed Orsiro BP-SES hybrid stent system. The trial included 30 patients with de novo lesions and demonstrated a low rate of in-stent restenosis with late lumen loss (LLL) of 0.05 ± 0.22 mm over the 9-month follow-up period. Though small, the BIOFLOW I study demonstrated that there was merit in the technology and, more importantly, that it was safe, with no myocardial infarction (MI) or stent thrombosis reported.

The follow-on BIOFLOW II [3] study was a prospective, international, randomized controlled, noninferiority trial that compared Orsiro BP-SES in a 2:1 allocation ratio with the benchmark stent of the day, Xience Prime (Abbott Vascular, Santa Clara, California). The trial also included a subgroup analysis, which allocated patients to repeat optical coherence tomography (OCT) for disease occurring in small vessels (<2.75 mm diameter) and diabetic patients ($n = 55$). The initial 9-month follow-up data published in 2015 showed non-inferior LLL compared with Xience (Orsiro 0.10 ± 0.32 mm vs. Xience 0.11 ± 0.29 mm, $p_{\text{noninferiority}} < 0.0001$). Similar to BIOFLOW I, no events of definite or probable stent thrombosis were reported in the Orsiro BP-SES group (298 patients, 332 lesions) [3]. The 5-year clinical follow-up data released in 2018 showed noninferior target lesion failure (TLF) (Orsiro 10.4% vs. Xience 12.7%, $p = 0.473$) with significantly lower mortality in vessels <2.75 mm (Orsiro 3.7% vs. 11.3%, $P = 0.022$), [10] suggesting that the 60- μ m thin Orsiro BP-SES profile may provide additional benefit in small-vessel subgroups. The recent publication of BIOFLOW VI [11] in China confirms the noninferior rates of LLL (Orsiro 0.05 mm \pm 0.21 mm vs. Xience Prime 0.07 mm \pm 0.2 mm) and TLF (Orsiro 2.3% vs. Xience 1.8%, $p = 0.7505$) when compared with Xience.

The BIOFLOW III [12,13] trial built on the findings of BIOFLOW II and aimed to gather “real world” data on TLF at 12 months using a prospective, open-label, registry-based design. Enrolling close to 1400 patients,

the registry reported a TLF rate of 5.1% over 12 months (95% CI 4.0% to 6.4%), with a rate of 10% (95% CI 8.4% to 12.0%) in the 5-year follow-up study [12]. Low rates of stent thrombosis similar to previous trials were observed (0.2% at 12 months, 95% CI 0.1% to 0.7%).

BIOFLOW IV [14] was designed as a noninferiority trial to allow regulatory approval in Japan. A total of 579 patients with de novo lesions were randomized in a 2:1 fashion to receive Orsiro BP-SES vs. the Xience permanent-polymer everolimus-eluting stent (PP-EES). Patients for planned intervention and patients presenting with MI in the previous 72 h were excluded. A total number of 385 patients (441 lesions) were enrolled. Twelve-month follow-up data reported a noninferior rate of target vessel failure (TVF) in de novo lesions compared with Xience (Orsiro 5.5% vs. Xience 7.5%, $p_{\text{noninferiority}} < 0.0001$). These data echoed the previous European BIOFLOW I, II, and III findings with similar rates of TLF and stent thrombosis, demonstrating applicability of those findings to the Japanese population and facilitating Japanese regulatory approval in January 2018 [15]. Five-year follow-up is ongoing, with results expected in 2024.

With regulatory approval in hand in both Europe and Japan, the BIOFLOW V [2] trial randomized patients to receive Orsiro BP-SES or Xience in patients presenting for elective or urgent percutaneous coronary intervention (PCI) at 90 hospitals across 13 countries, including the United States for the first time. The primary endpoint of the trial was TLF at 12 months. Of patients included in the study, 51% presented with an acute coronary syndrome. Although BIOFLOW V was powered as a non-inferiority trial, over the initial 12-month period, the incidence of TLF was significantly lower in the Orsiro BP-SES group than in the Xience DP-EES group (Orsiro 6.2% vs. Xience 9.6%, $p = 0.0399$) and was likely driven by a lower rate of target vessel MI in the Orsiro BP-SES group (Orsiro 5% vs. Xience 8%, $p = 0.0155$) [16].

In order to improve endpoint detection and statistical significance, the authors of BIOFLOW V incorporated the results of BIOFLOW II and IV and used a Bayesian analysis [17]. This pooled analysis showed a TLF rate of 3% in the Orsiro BP-SES group with a posterior probability for noninferiority of 100% and a posterior probability of superiority of 97% [2]. A subsequent landmark analysis of the BIOFLOW V data confirmed that the statistically significant differences in TLF and target vessel MI reported at one year persisted through two-year follow-up [18]. Shortly after BIOFLOW V's publication, the Food and Drug Administration granted approval for use of the device in the United States [19].

2.2. The BIOSCIENCE, BIO-RESORT, BIONYX, and SORT OUT VII and IX trials

In addition to the data provided by the BIOFLOW studies, several large-scale investigator-initiated trials demonstrating the safety and efficacy of the Orsiro BP-SES have also emerged. Fig. 3 shows rates of TLF at 12 months from selected prospective randomized controlled trials comparing the Orsiro BP-SES with the Xience PP-EES.

The BIOSCIENCE [1] trial, published in 2014, was similar to BIOFLOW II in design and was a large scale, 1:1 randomized noninferiority trial. Unlike BIOFLOW II, however, patients with acute coronary syndrome in the preceding 72 h were included. Patient enrollment was completed in 12 months, with 2119 patients and 3139 lesions enrolled. The primary endpoint of TLF at 12 months showed that Orsiro BP-SES was non-inferior to Xience DP-EES (Orsiro 6.5% vs. Xience 6.6%, $p_{\text{noninferiority}} < 0.0004$), similar to the findings of BIOFLOW II. A novel finding, however, was the subgroup analysis in patients treated with ST-elevation myocardial infarction (STEMI). In the prespecified subgroup of acute coronary syndromes with STEMI, the Orsiro BP-SES was associated with a lower risk of TLF than the Xience PP-EES at 12 months (Orsiro 3.3% vs. Xience 8.7%, $p = 0.024$) [1,20]. The BIOSCIENCE 5-year follow-up data [20] failed to show any difference in the TLF rate at 5 years between the two groups, suggesting that the reduction seen at 12 months may be related to the BP-SES's ultrathin struts, and that after complete











Fig. 1. A) The open-cell double-helix of the ultrathin-strut Orsiro biodegradable-polymer sirolimus-eluting stent is designed for flexibility and deliverability. B) Cross-section of the ultrathin biodegradable-polymer strut showing the innermost cobalt chromium alloy PRO-Kinetic energy™ stent, the middle proBIO™ silicon carbide layer that seals the metal alloy surface and reduces the allergenic interaction of metal ions with the vessel wall and the outer BIOlute™ layer is composed of a bioabsorbable poly-L-lactic acid (PLLA) polymer containing the antiproliferative sirolimus agent. Images provided by and used with the permission of Biotronik (BIOTRONIK, Buelach, Switzerland).

breakdown of the biodegradable polymer, the risk of very late TLF is similar to that observed with permanent-polymer stents.

One question that remained after the publication of the early BIOFLOW trials and the BIOSCIENCE trial was whether the low rates of TLF observed with the Orsiro BP-SES were due to the thin strut design and the absence of a permanent-polymer stent within the vessel wall, reducing vascular inflammation and the development of

neoatherosclerosis. The BIO-RESORT [21] study aimed to answer this question using a large-scale, multicenter three-arm randomized control trial. The 3514 patients who were enrolled included all-comers and acute coronary syndromes with minimal exclusion factors and randomized patients to receive one of three stents. The very-thin-strut Orsiro BP-SES (60 μm), the very-thin-strut Synergy biodegradable-polymer everolimus-eluting stent (74 μm , BP-EES) or the thin-strut Resolute

Abbott/Boston Xience/Promus	Medtronic Resolute Onyx	Biosensors BioMatrix	Terumo Nobori	Terumo Ultimaster	Boston Synergy	Biotronik Orsiro
STENT PLATFORM						
						
CoCr/PtCr	CoNi	316L SS	316L SS	CoCr	PtCr	CoCr
STRUT THICKNESS						
81 μm	81 μm	120 μm	120 μm	80 μm	74 μm	60 μm
POLYMER COATING						
DURABLE POLYMER		BIOABSORBABLE POLYMER				
PVDF-HFP	BioLinx	PLA	PLA	PDLA-PCL copolymer	PLGA	PLLA
POLYMER DEGRADATION						
n/a	n/a	6-9 mo	6-9 mo	3-4 mo	4 mo	> 12 mo
POLYMER DISTRIBUTION						
Circumferential	Circumferential	Abluminal	Abluminal	Abluminal	Abluminal	Circumferential
7-8 $\mu\text{m}/\text{side}$	6 $\mu\text{m}/\text{side}$	10 μm	20 μm	15 μm	4 μm	4-7 $\mu\text{m}/\text{side}$
DRUG TYPE/DOSE						
Everolimus	Zotarolimus	Biolimus A9	Biolimus A9	Sirolimus	Everolimus	Sirolimus
1.0 $\mu\text{g}/\text{mm}^2$	1.6 $\mu\text{g}/\text{mm}^2$	15.6 $\mu\text{g}/\text{mm}$	15.6 $\mu\text{g}/\text{mm}$	3.9 $\mu\text{g}/\text{mm}$	113 $\mu\text{g} / 20 \text{ mm}$	1.4 $\mu\text{g}/\text{mm}^2$

Adapted from Iglesias et al., Orsiro cobalt-chromium sirolimus-eluting stent: present and future perspectives, DOI: 10.1080/17434440.2017.1378091

Fig. 2. Figure displaying the comparative strut thicknesses, polymer coatings, polymer degradation periods, distribution of the polymer, the eluted drug and drug load. Images provided by and used with the permission of Biotronik (BIOTRONIK, Buelach, Switzerland), from: Iglesias et al., Orsiro cobalt-chromium sirolimus-eluting stent: present and future perspectives, DOI: <https://doi.org/10.1080/17434440.2017.1378091>

Table 1

Summary table of studies describing study design highlights, number of patients enrolled, number of lesions treated, primary endpoints, and key findings of large scale, prospective randomized controlled trials and registries for the Orsiro bioresorbable-polymer sirolimus-eluting stent (BP-SES). LLL = Late lumen loss, TLF = Target lesion failure (composite of cardiovascular death, target vessel-related myocardial infarction, or ischemia-driven target lesion revascularization), TVR = Target vessel failure (composite of cardiac death, target vessel-related myocardial infarction, or clinically indicated target vessel revascularization), MACE = major adverse cardiovascular events, TLR = target vessel revascularization, MI = myocardial infarction, OCT = optical coherence tomography.

Study name	Study design highlights	Total patients	Patients in BP-SES arm	Lesions in BP-SES arm	Primary endpoint	Key findings
BIOFLOW I	- First in man - Multi-center - Single arm registry - Single de novo lesions	30	30	30	LLL at 9 months	- LLL 0.05 ± 0.22 mm at 9 months. - 10% MACE (1 cardiac death, 2 non-cardiac deaths). - 2 patients with TLR. - No MI or stent thrombosis.
BIOFLOW II - 12 month	- Orsiro vs. Xience Prime - ≤ 2 De-novo lesions	452	298	332	LLL at 9 months	- Non-inferior LLL at 9 months - Orsiro 0.10 ± 0.32 mm vs. Xience 0.11 ± 0.29 mm (p -value for non-inferiority <0.0001). - Clinical endpoints comparable to Xience over 12-month follow up.
BIOFLOW II - 5 year					TLF at 60 months	- Non-inferior TLF at 5 years - Orsiro 10.4% vs. Xience 12.7% ($p = 0.473$). - No stent thrombosis in Orsiro group. - Numerically more TLR in diabetes treated with Orsiro, but not statistically significant (13.5% vs. 4.5%, $p = 0.138$). - Orsiro showed significantly lower 5-year mortality in small vessel disease (<2.75 mm) - (3.7% vs. 11.3%, $p = 0.022$). - TLF rate of 5.1% at 12 months (95% CI 4.0–6.4%). - Cardiac death 1.3% (95% CI 0.9–2.1%). - Definite stent thrombosis 0.2% (95% CI 0.1–0.7%).
BIOFLOW III	- Single arm registry - All-comers	1356	1356	1738	TLF at 12 months	- TLF rate of 10% at 5 years (95% CI 8.4–12%). - TLF, cardiac death and MI rates higher in diabetes compared with non-diabetes, but TLF still low at 14%
BIOFLOW III - 5 year					TLF at 60 months	- Non-inferior TVF at 12 months - Orsiro 5.5% vs. Xience 7.5% (p -value for non-inferiority <0.0001). - 5-year follow up is ongoing.
BIOFLOW IV	- Orsiro vs. Xience Prime (2:1) - ≤ 2 De-novo lesions - ACS excluded	575	385	441	TVF at 12 months	- Superior rate of TLF compared to Xience* - Orsiro 6.2% vs. 9.6% in Xience group at 12 months ($p = 0.0399$). - Lower cardiac death or MI in Orsiro group at 24 months (5% vs. 9%, $p = 0.072$).
BIOFLOW V	- Orsiro vs. Xience (2:1) - STEMI, bypass grafts, ISR, CTO's excluded	1334	884	–	TLF at 12 months	- Non-inferior LLL at 9 months - Orsiro 0.05 mm ± 0.21 mm vs. Xience Prime 0.07 mm ± 0.2 mm (p -value non-inferiority <0.0001). - TLF similar - Orsiro 2.3% vs. Xience 1.8% ($p = 0.7505$) - Final publication pending
BIOFLOW VI	- Orsiro vs. Xience Prime (1:1) - ACS, grafts, CTO's, bifurcation & calcified lesions excluded	440	223	257	LLL at 9 months	- Non-inferior to DCB in ISR - LLL in Orsiro 0.2 mm vs. 0.03 mm in DCB group ($p = 0.39$).
BIOLUX RCT	- Pantera Lux drug coated balloon (DCB) with Orsiro in ISR (2:1)	229	72	80	LLL at 6 months	- Non-inferior TLF at 12 months - Orsiro 6.5% vs. Xience 6.6% (p -value for non-inferiority <0.0004)
BIOSCIENCE	- Orsiro to Xience Prime (1:1) - All-comers	2119	1063	1594	TLF at 12 months	- Non-inferior TLF at 5 years - Orsiro 20.2% vs. Xience 18.8%. In subgroup analysis Orsiro demonstrated lower TLF in STEMI patients at 5 years compared to Xience (RR 0.38, 0.16–0.91)
BIOSCIENCE - 5 year					TLF at 60 months	- Non-inferior TVF at 12 months - (Orsiro 4.7% vs. Resolute Integrity 5.4%, p non-inferiority <0.0001) - Orsiro had numerically lowest rates of TVF, TLF, TLR, cardiac death, stent thrombosis and MACE - Landmark analysis suggests Orsiro may reduce risk of TLR at 1 year
BIO-RESORT	- Orsiro, Synergy or Resolute Integrity (1:1:1) - All-comers	3514	1169	1551	TVF at 12 months	- Non-inferior TLF at 12 months - Orsiro 3.8 vs. Nobori 4.6% (0.12–0.92; $P = 0.034$). - Orsiro demonstrated reduced risk of definite stent thrombosis at 12 months (0.4% vs. 1.2%, $p = 0.03$). - Non-inferior TVF at 12 months - Orsiro 4.5% vs. Resolute Onyx 4.7% (p -value for non-inferiority <0.0005). - Definite or probable stent thrombosis low in both groups - Orsiro 0.7% vs. Resolute Onyx 0.1%. - BioFreedom did not meet criteria for non-inferior TLF compared with Orsiro (Orsiro 4.0% vs. BioFreedom 5.3%, p -value non-inferiority <0.123). - Orsiro had significantly less TLR at 12 months (Orsiro 1.3% vs. BioFreedom 3.5%, $p \leq 0.0001$). - Ongoing
SORT OUT VII	- Orsiro vs. Nobori (1:1) - All-Comers	2525	1261	1590	TLF at 12 months	- Superior rates of TLF at 12 months - Orsiro 4% vs. Xience 6% (rate ratio 0.59, posterior probability of superiority 0.986)
BIONYX	- Orsiro vs. Resolute Onyx (1:1) - All-comers	2488	1245	1593	TVF at 12 months	
SORT OUT IX	- Orsiro vs. BioFreedom (1:1) - All-comers	3150	1579		TLF at 60 months	
SORT OUT X	- Orsiro vs. CD4 Combo stent (1:1)	3148	c.1570	–	TLF at 12 months	
BIOSTEMI	- Superiority trial - Orsiro vs. Xience	1300	649	–	TLF ^U at 12 months	

(continued on next page)

Table 1 (continued)

Study name	Study design highlights	Total patients	Patients in BP-SES arm	Lesions in BP-SES arm	Primary endpoint	Key findings
ORIENT	(1:1) - STEMI patients - Orsiro vs. Resolute Integrity (2:1) - All-comers	372	250	345	LLL at 9 months	- Non-inferior LLL at 9 months - (Orsiro 0.10 ± 0.35 mm, Resolute Integrity 0.16 ± 0.39 mm, p for non-inferiority <0.001)
HATTRICK-OCT	- Orsiro vs. Endeavour Resolute (1:1) - ACS patients with de novo lesions in LAD	44	23	23	Strut coverage by OCT, vasodilator response at 3 months	- Proportion of uncovered struts by OCT less in Orsiro (Orsiro 3.9% vs. Endeavour Resolute 8.9%) at 3 months.
PRISON IV	- Orsiro vs. Xience CTO's (1:1)	330	165	165	LLL at 9 months	- Orsiro failed to meet criteria for non-inferior LLL at 9 months in recanalized CTO's (Orsiro 0.13 ± 0.63 mm vs. Xience 0.02 ± 0.47 mm, p non-inferiority = 0.11), sub group analysis suggested difference most pronounced in Orsiro stents <3.0 mm.
		23,406	10,652	9739		

Integrity durable-polymer zotarolimus-eluting stent (91 μ m, DP-ZES). Over 12 months, the Orsiro BP-SES proved to be noninferior to the durable-polymer Resolute Integrity stents and had similar rates of TVF as the Synergy BP-EES [21,22].

Follow-up data at 2 years showed that TVF was not statistically different between the three stents; however, when TLF was examined using a landmark analysis between 12 and 24 months, Orsiro BP-SES fared better than Resolute, with a target vessel revascularization (TVR) rate of 0.6% for Orsiro compared with 1.5% for Resolute, suggesting that Orsiro BP-SES reduced the risk of revascularization after 1 year of follow-up [22].

Based on the rates of TLF seen in the BIOFLOW V [2] and BIO-RESORT [21,22] studies, the Orsiro BP-SES was selected as the comparator for the first randomized controlled trial of the newly designed Resolute Onyx stent (Medtronic, Santa Rosa, California, USA) in the BIONYX trial [23]. Resolute Onyx was developed as a thin-strut durable-polymer zotarolimus-eluting stent system. The trial included all comers, randomizing 2516 patients in a 1:1 fashion. Both stents showed excellent safety and efficacy, with low rates of definite or probable stent thrombosis in both groups (Orsiro 0.7% vs. Onyx 0.1%, hazard ratio 0.11 [95% CI 0.01–0.87]; $p = 0.0112$). The primary endpoint of TLF at 12 months was similar in both the Orsiro BP-SES and Onyx groups [23]. Funding has been secured for up to three years of follow-up, with two-year data expected in the coming months.

Whereas the BIOFLOW, BIO-RESORT, BIOSCIENCE, and BIONYX trials compared Orsiro BP-SES to durable-polymer systems, the SORT OUT VII trial [24] was the first head-to-head trial against a biodegradable-polymer system. The registry-based, multicenter, noninferiority trial randomized patients to receive the Orsiro BP-SES or the Nobori (stainless-steel, biolimus-eluting, 120- μ m strut) stent. The trial enrolled 2525 patients, with TLF at 12 months being the primary endpoint. Orsiro BP-SES was shown to be noninferior to Nobori at 12 months with respect to TLF (Orsiro 3.8% vs. Nobori 4.6%, 95% CI 0.12–0.92; $p_{\text{noninferiority}} = 0.034$). The rates of reported TLF were most likely lower than reported in the large BIOSCIENCE trial (6.5%), as the SORT OUT (3.8%) trial excluded procedure-related MI because of the registry-based design. In addition, the Orsiro BP-SES demonstrated a statistically lower rate of definite or probable stent thrombosis at 12 months (Orsiro 0.4% vs. Nobori 1.2%, $p = 0.03$) [24].

Further studies against polymer-free drug-coated stents, such as the SORT OUT IX [25] trial, showed that the biolimus-based BioFreedom stent failed to meet the criteria for noninferiority of TLF at 12 months compared with Orsiro BP-SES, likely driven by a lower rate of target lesion revascularization (TLR) in the Orsiro BP-SES arm of the study [26]. Recent meta-analyses have suggested that the improved rates of TLF may be related to an increased risk of stent thrombosis in thick-strut platforms such as 120- μ m Nobori and 120- μ m BioFreedom stents compared with the 60- μ m Orsiro BP-SES system [27–29].

Data provided by the Swedish Coronary Angiography and Angioplasty Registry (SCAAR) corroborate the results seen in randomized trials. An analysis over the 6-year period from October 2011 to June 2017 including 74,131 patients (4561 in the Orsiro group compared to 69,570 who received other frequently used DES) showed significantly lower rates of TLR in the Orsiro group (Orsiro 1.6% vs. other DES 2.3%, adjusted HR 0.75, 95% CI 0.60–0.94; $p = 0.012$) [30].

2.3. Evidence in high-risk subgroups

PCI in high-risk subgroups can present several challenging patient and lesion characteristics and can result in increased rates of in-stent restenosis, TVF, and stent thrombosis. The following section will address the performance of the Orsiro BP-SES in these high-risk subgroups.

2.3.1. ST-elevation myocardial infarction

Patients presenting with STEMI pose several challenges to PCI. Malapposition of stent struts may result from undersizing and subsequent thrombus resolution, setting the stage for an increased risk of early stent thrombosis [31]. Although the development of DES allowed for substantial improvements in early stent thrombosis, concerns remained about delayed vascular healing due to hypersensitivity and inflammation induced by permanent-polymer stents [32,33]. The thin-strut profile coupled with the inert silicon carbide coating of the Orsiro BP-SES has shown promise in reducing rates of TLF, based on subgroup analyses in larger recent trials [2,20,31,34].

The BIOSCIENCE trial's 12- and 24-month [1,22] results demonstrated improved outcomes in STEMI patients. Twelve-month rates of TLF were reported at 3.3% vs. 8.7% in the BP-SES and DP-EES groups (risk ratio [RR] 0.38, $p = 0.024$) with 24-month data showing rates of 5.4% vs 10.8% (RR 0.48, $p = 0.043$) [32]. Most notably, the lower observed rates of TLF appeared to be driven by decreased rates of cardiac death or MI in the BP-SES group (RR 0.46, 95% CI 0.21–1.02, $p = 0.05$) [32]. BIOFLOW V data showed a similar reduction in TLF with a reported hazard ratio of 0.50 (95% CI 0.28–0.89) compared to DP-EES [2].

The combined results of the BIOFLOW V and BIOSCIENCE noninferiority trials and the demonstrated improved outcomes led to the design of the BIOSTEMI trial [35]. The trial enrolled 1300 patients presenting with STEMI with 1:1 randomization to the BP-SES or DP-EES treatment group. Unlike previous trials powered for noninferiority, BIOSTEMI demonstrated, for the first time, clear superiority with respect to rates of TLF at 12 months in patients presenting with STEMI (Orsiro 4% vs. Xience 6%, rate ratio 0.59, posterior probability of superiority 0.986) [36]. The difference in the rates of TLF was driven primarily by fewer cases of clinically indicated TLR in the Orsiro group.




Fig. 3. Rates of target lesion failure at 12 months from selected prospective randomized controlled trials comparing the Orsiro biodegradable-polymer sirolimus eluting stent (BP-SES) to the Xience permanent-polymer everolimus-eluting stent (PP-EES).

2.3.2. Chronic total occlusions

Chronic total occlusions (CTOs) present one of the toughest environments for stent performance with respect to in-stent thrombosis and late or very late in-stent restenosis [37,38]. When compared to de novo lesions, CTOs are typically longer and more calcified, often require multiple stents with overlapping segments, and may undergo late vessel remodeling as blood flow increases, leading to vasodilation and stent-strut malapposition, setting the stage for thrombosis or in-stent restenosis [39,40]. The PRISON IV [40] trial set out to evaluate the Orsiro BP-SES performance in this challenging population. Enrolling 330 patients, the trial randomized patients to receive Orsiro BP-SES or Xience in a 1:1 allocation. Of patients enrolled in both arms of the trial, 99% underwent repeat angiography at 9 months and clinical follow-up at 12 months. PRISON IV failed to meet the primary endpoint of noninferior LLL at 9 months of Orsiro BP-SES compared with Xience and showed a significantly higher rate of binary restenosis in the Orsiro BP-SES group (Orsiro 8.0% vs. Xience 2.1%, $p = 0.028$).

The recent publication of the 3-year follow up data of the PRISON IV trial [41] gives the best long-term follow up of the Orsiro BP-SES's performance in CTOs. Over the three years of follow-up, the cumulative incidence of major adverse clinical events was higher in the Orsiro BP-SES arm than in the Xience EES arm (9.7% vs. 4.2%). Similar trends were also seen with respect to TLR over the study period favoring Xience with TLR rates reported as 9.7% in the Orsiro arm versus 3.4% in the Xience arm. A landmark analysis demonstrated that after one year, rates of TLR overall remained low but were higher in the Orsiro arm than in the Xience arm (2.1% vs. 0.6%, p -value reported as "NS") [41].

In a smaller, non-randomized study, Markovic et al. [38] examined LLL at 9 months in patients receiving the Resolute Integrity zotarolimus-eluting stent versus the Orsiro BP-SES. Rates of LLL were significantly less in the Orsiro BP-SES arm (Orsiro 0.24 ± 0.53 mm vs. Resolute Integrity 0.59 ± 0.72 mm, $p = 0.01$), but this did not translate into a difference in clinical events up to 24 months of follow-up.

Although PRISON IV failed to show noninferiority with respect to LLL, the study was underpowered to detect clinical endpoints. The rates of TVR were higher in PRISON IV compared with BIOSCIENCE (9.2% vs. 4%), but this was probably due to the routine angiographic follow-up at 9 months in the PRISON IV trial as opposed to reported clinical events in the BIOSCIENCE trial. Subgroup analyses have demonstrated that the majority of the LLL seen within the PRISON IV study is due to stents with <3 mm in diameter (60- μ m strut thickness) vs. stents >3 mm in diameter (80- μ m strut thickness) [42]. An additional OCT subgroup analysis demonstrated favorable strut coverage in the Orsiro arm [43]. Despite these findings, results from CTO patients in BIOFLOW III

[12,13] and SORT OUT VII [44] suggest that rates of TLF and TLR remain low in CTO subgroups.

2.3.3. Diabetes

A recent meta-analysis pooling diabetic patients from the BIOFLOW II, IV, and V trials examined TLF at 1 year. The analysis included 494 diabetic patients treated with Orsiro BP-SES and 263 patients treated with Xience across the three BIOFLOW trials. Rates of TLF were similar at one year: 6.3% in the Orsiro BP-SES group and 8.7% in the Xience group (HR 0.82, 95% CI 0.047–1.43, $p = 0.493$) [16]. These findings were confirmed in a subgroup analysis of the SORT OUT VII trial [45]. Rates of TLF among patients treated with or without insulin were also found to be similar [45].

2.3.4. Small-vessel disease

Small-vessel disease, generally defined as vessels <2.75 mm in diameter, has been associated with a greater risk of TLR due to increased rates of in-stent restenosis [46,47]. Data from BIOFLOW II [10] and a subgroup analysis from BIO-RESORT [48] suggest that the ultrathin Orsiro BP-SES may provide additional benefit in small-vessel subgroups.

BIOFLOW II contained 259 patients within the small-vessel subgroup (≤ 2.75 mm) and demonstrated 5-year TLF rates of 11.1% in the BP-SES arm vs. 15.5% in the DP-EES arm ($p = 0.303$), a difference that was mainly due to a lower rate of death in the BP-SES arm (3.7% vs. 11.3%, $p = 0.039$) [10]. A subsequent analysis of the BIO-RESORT trial examined the 3-year outcomes in 1506 of the 3514 total trial participants. A multivariate analysis later showed that BP-SES was independently associated with lower rates of TLR at 3 years in comparison with the everolimus and zotarolimus treatment groups (adjusted HR 0.42, 95% CI 0.20–0.85, $p = 0.02$). No differences were reported in the rates of cardiac death, target vessel MI, or stent thrombosis at 3 years [48].

Results from the small-vessel cohort of the BIOSCIENCE trial [47] released in August 2019 differed somewhat from the results of the BIOFLOW II and BIO-RESORT small-vessel subgroups. The ≤ 2.75 -mm small-vessel cutoff used in BIOFLOW II and BIO-RESORT, versus the definition of ≤ 3 mm or > 3 mm to define small- versus large-vessel disease used in BIOSCIENCE, potentially accounts for the discrepant results. Five-year rates of TLF were not found to be significantly different between the BP-SES and DP-EES groups (22.3% vs. 18.3%, respectively; rate ratio 1.22, 95% CI 0.94–1.58, $p = 0.13$) [47].

2.3.5. In-stent restenosis

Despite recent advances in stent technology, neoatherosclerosis due to inflammation, drug toxicity, and delayed vascular healing can result in revascularization rates of up to 10% [4,49]. The inert silicon carbide proBIO layer has the potential to improve biocompatibility at the vessel interface and to reduce the proliferation of tissue, which leads to in-stent restenosis. However, to date, there are limited available randomized data on the efficacy of BP-SES versus existing second-generation DES in treating in-stent restenosis, with only 10 patients within the all-comer BIOSCIENCE trial and 30 within BIO-RESORT.

Current European guidelines advocate for the use of drug-coated balloons (DCBs) or DES in the treatment of in-stent restenosis, but there are limited data on the efficacy of latest-generation stents [50]. The BIOLUX trial [49] aimed to address this gap, randomizing 229 patients with in-stent restenosis in a 2:1 fashion to receive treatment with a butyryl trihexyl citrate (BTHC)-based paclitaxel DCB versus BP-SES. In-stent LLL at 6 months and TLF at 12 months were compared between both groups. Over the study period, DCB proved noninferior to BP-SES, with LLL in DCB arm reported at 0.03 ± 0.40 mm and 0.20 ± 0.70 mm in the BP-SES arm ($p = 0.40$). Rates of TLF at 12 and 18 months were also similar, demonstrating that both options seem feasible and effective but that further studies are needed in this challenging cohort of patients [49].

2.4. Economic impact

In addition to the improved clinical outcomes demonstrated in the BIOFLOW V trial, recent economic modeling has also suggested possible economic benefits [51]. Markov economic modeling has shown that the reduction in rates of TLR and target vessel MI resulted in a cumulative reduction of \$2415 per patient over a 48-month period compared to Xience [52]. With rising healthcare costs, these savings are an important consideration, not only for individual patients and insurance payers, but also for healthcare policymakers [51]. While caution should be exercised in interpreting this single sub-analysis, healthcare spending in the US alone is projected to reach 20% of gross domestic product by 2025 [53], so the delivery of safe, clinically effective and potentially economically viable treatments is an important consideration for any healthcare system.

3. Limitations

The Orsiro BP-SES does have some notable limitations. The presence of the ultrathin struts may make the stent a less favorable choice in situations requiring high radial force to prevent vessels recoil, such as CTOs, as described in a previous section. Additionally, data on the performance of the Orsiro BP-SES in heavily calcified lesions is limited, mainly arising from small subgroups of “all-comers” trials, making any definite statements about the device's performance in this space difficult.

On a practical level, the ultrathin-stent struts can be difficult to visualize using fluoroscopy, so careful confirmation of the stent position should be performed before deployment.

4. Future directions

Several ongoing trials are continuing to investigate the efficacy and effectiveness of the Orsiro BP-SES in comparison to other contemporary stents. The enhanced biocompatibility and ultrathin-strut design may also have a role to play in high-bleeding-risk populations by allowing for shortened dual antiplatelet therapy duration, which will be examined in the planned BIOFLOW-DAPT study. With the favorable outcomes of the Orsiro BP-SES, ongoing developments in stent technology, such as the incorporation of CD34+ antibodies into bioabsorbable-polymer stents (COMBO™, OrbusNeich, Hoevalaken, the Netherlands), aimed

at encouraging endothelial progenitor cell growth and vascular healing, are now being compared to the Orsiro BP-SES [54].

5. Conclusions

Overall, the Orsiro BP-SES has demonstrated its safety and noninferiority in numerous randomized clinical trials [1,3,14,20,48]. The combination of the ultrathin-strut design and inert proBIO silicon carbide layer allowing enhanced biocompatibility has demonstrated low rates of TLR against comparators and low rates of stent thrombosis. In specific high-risk subgroups, such as patients with STEMI and CTOs, notable exceptions have emerged, with lower rates of cardiac death and MI in STEMI patients [36], yet higher rates of TLR and major adverse cardiac events suggested in patients with CTOs [41]. Additional appropriately powered randomized controlled studies are needed to examine the role of the Orsiro BP-SES in other high-risk subgroups, including high-bleeding-risk patients, CTOs, and in-stent restenosis.

Funding

This research did not receive any specific grant from funding agencies in the public, commercial, or not-for-profit sectors.

Declaration of competing interest

Ron Waksman – Advisory Board: Amgen, Boston Scientific, Cardioset, Cardiovascular Systems Inc., Medtronic, Philips, Pi-Cardia Ltd.; Consultant: Amgen, Biotronik, Boston Scientific, Cardioset, Cardiovascular Systems Inc., Medtronic, Philips, Pi-Cardia Ltd.; Grant Support: AstraZeneca, Biotronik, Boston Scientific, Chiesi; Speakers Bureau: AstraZeneca, Chiesi; Investor: MedAlliance.

All other authors: None.

References

- [1] Pilgrim T, Heg D, Roffi M, Tuller D, Muller O, Vuilliamet A, et al. Ultrathin strut biodegradable polymer sirolimus-eluting stent versus durable polymer everolimus-eluting stent for percutaneous coronary revascularisation (BIOSCIENCE): a randomised, single-blind, non-inferiority trial. *Lancet* 2014;384(9960):2111–22.
- [2] Kandzari DE, Mauri L, Koolen JJ, Massaro JM, Doros G, Garcia-Garcia HM, et al. Ultrathin, bioresorbable polymer sirolimus-eluting stents versus thin, durable polymer everolimus-eluting stents in patients undergoing coronary revascularisation (BIOFLOW V): a randomised trial. *Lancet* 2017;390(10105):1843–52.
- [3] Windecker S, Haude M, Neumann FJ, Stangl K, Witzentichler B, Slagboom T, et al. Comparison of a novel biodegradable polymer sirolimus-eluting stent with a durable polymer everolimus-eluting stent: results of the randomized BIOFLOW-II trial. *Circ Cardiovasc Interv* 2015;8(2):e001441.
- [4] Iglesias JF, Muller O, Zuffi A, Eeckhout E. Performance of the Orsiro hybrid drug-eluting stent in high-risk subgroups. *Minerva Cardioangiol* 2016;64(1):55–73.
- [5] Tittelbach M, Diener T. Orsiro - the first hybrid drug-eluting stent, opening up a new class of drug-eluting stents for superior patient outcomes. *Interventional Cardiology Review* 2011;6(2):142–4.
- [6] Hamon M, Niculescu R, Deleanu D, Dorobantu M, Weissman NJ, Waksman R. Clinical and angiographic experience with a third-generation drug-eluting Orsiro stent in the treatment of single de novo coronary artery lesions (BIOFLOW-I): a prospective, first-in-man study. *EuroIntervention* 2013;8(9):1006–11.
- [7] Finn AV, Kolodgie FD, Harnek J, Guerrero LJ, Acampado E, Tefera K, et al. Differential response of delayed healing and persistent inflammation at sites of overlapping sirolimus- or paclitaxel-eluting stents. *Circulation* 2005;112(2):270–8.
- [8] Joner M, Finn AV, Farb A, Mont EK, Kolodgie FD, Ladich E, et al. Pathology of drug-eluting stents in humans: delayed healing and late thrombotic risk. *J Am Coll Cardiol* 2006;48(1):193–202.
- [9] Onuma Y, Serruys P, den Heijer P, Joesoef KS, Duckers H, Regar E, et al. MAHOROBA, first-in-man study: 6-month results of a biodegradable polymer sustained release tacrolimus-eluting stent in de novo coronary stenoses. *Eur Heart J* 2009;30(12):1477–85.
- [10] Lefevre T, Haude M, Neumann FJ, Stangl K, Skurk C, Slagboom T, et al. Comparison of a novel biodegradable polymer sirolimus-eluting stent with a durable polymer everolimus-eluting stent: 5-year outcomes of the randomized BIOFLOW-II trial. *JACC Cardiovasc Interv* 2018;11(10):995–1002.
- [11] Yang Y, editor. BIOFLOW VI -Oral Presentation CIT 2018, 2018. ; Mar 23 2018 Suzhou, China. Available: <http://www.citmd.com/CIT/2018/generic/web/viewer.php?strencode=f18Tnhw1h8vWylCPGV9Qj5ddEl8dSIDdGMuTCh2Bh5/TfDMenchPw2mMkYEZgBSNGd3RWtGWWQs=&spid=560>.

- [12] Waltenberger J, Brachmann J, van der Heyden J, Richardt G, Frobert O, Seige M, et al. Five-year results of the BIOFLOW-III registry: real-world experience with a biodegradable polymer sirolimus-eluting stent. *Cardiovasc Revasc Med* 2020;21(1):63–9.
- [13] Waltenberger J, Brachmann J, van der Heyden J, Richardt G, Frobert O, Seige M, et al. Real-world experience with a novel biodegradable polymer sirolimus-eluting stent: twelve-month results of the BIOFLOW-III registry. *EuroIntervention* 2016;11(10):1106–10.
- [14] Saito S, Toelg R, Witzensbichler B, Haude M, Masotti M, Salmeron R, et al. BIOFLOW-IV, a randomised, intercontinental, multicentre study to assess the safety and effectiveness of the Orsio sirolimus-eluting stent in the treatment of subjects with de novo coronary artery lesions: primary outcome target vessel failure at 12 months. *EuroIntervention* 2019;15(11):e1006–e13.
- [15] Notarangelo FM, Maglietta G, Bevilacqua P, Cereda M, Merlini PA, Villani GQ, et al. Pharmacogenomic approach to selecting antiplatelet therapy in patients with acute coronary syndromes: the PHARMICLO trial. *J Am Coll Cardiol* 2018;71(17):1869–77.
- [16] Waksman R, Shlofmitz E, Windecker S, Koolen JJ, Saito S, Kandzari D, et al. Efficacy and safety of ultrathin, Bioreabsorbable-polymer Sirolimus-eluting stents versus thin, durable-polymer Everolimus-eluting stents for coronary revascularization of patients with diabetes mellitus. *Am J Cardiol* 2019;124(7):1020–6.
- [17] Doros G, Massaro JM, Kandzari DE, Waksman R, Koolen JJ, Cutlip DE, et al. Rationale of a novel study design for the BIOFLOW V study, a prospective, randomized multicenter study to assess the safety and efficacy of the Orsio sirolimus-eluting coronary stent system using a Bayesian approach. *Am Heart J* 2017;193:35–45.
- [18] Kandzari DE, Koolen JJ, Doros G, Massaro JJ, Garcia-Garcia HM, Bennett J, et al. Ultrathin bioreabsorbable polymer sirolimus-eluting stents versus thin durable polymer Everolimus-eluting stents. *J Am Coll Cardiol* 2018;72(25):3287–97.
- [19] United States Food and Drug Administration. Premarket Approval (PMA) - ORSIRO Sirolimus Eluting Coronary Stent System 2019. Available from <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpma/pma.cfm?id=P170030>.
- [20] Pilgrim T, Piccolo R, Heg D, Roffi M, Tuller D, Muller O, et al. Ultrathin-strut, biodegradable-polymer, sirolimus-eluting stents versus thin-strut, durable-polymer, everolimus-eluting stents for percutaneous coronary revascularisation: 5-year outcomes of the BIOSCIENCE randomised trial. *Lancet* 2018;392(10149):737–46.
- [21] von Birgelen C, Kok MM, van der Heijden LC, Danse PW, Schotborgh CE, Scholte M, et al. Very thin strut biodegradable polymer everolimus-eluting and sirolimus-eluting stents versus durable polymer zotarolimus-eluting stents in allcomers with coronary artery disease (BIO-RESORT): a three-arm, randomised, non-inferiority trial. *Lancet* 2016;388(10060):2607–17.
- [22] Kok MM, Zocca P, Buiten RA, Danse PW, Schotborgh CE, Scholte M, et al. Two-year clinical outcome of all-comers treated with three highly dissimilar contemporary coronary drug-eluting stents in the randomised BIO-RESORT trial. *EuroIntervention* 2018;14(8):915–23.
- [23] von Birgelen C, Zocca P, Buiten RA, Jessurun GAJ, Schotborgh CE, Roguin A, et al. Thin composite wire strut, durable polymer-coated (resolute onyx) versus ultrathin cobalt-chromium strut, bioreabsorbable polymer-coated (Orsio) drug-eluting stents in allcomers with coronary artery disease (BIONYX): an international, single-blind, randomised non-inferiority trial. *Lancet* 2018;392(10154):1235–45.
- [24] Jensen LO, Thayssen P, Maeng M, Ravkilde J, Kruse LR, Raugaard B, et al. Randomized comparison of a biodegradable polymer ultrathin strut sirolimus-eluting stent with a biodegradable polymer biolimus-eluting stent in patients treated with percutaneous coronary intervention: the SORT OUT VII trial. *Circ Cardiovasc Interv* 2016;9(7).
- [25] Jensen LO, Maeng M, Raugaard B, Engstrom T, Hansen HS, Jensen SE, et al. Comparison of the polymer-free biolimus-coated BioFreedom stent with the thin-strut biodegradable polymer sirolimus-eluting Orsio stent in an all-comers population treated with percutaneous coronary intervention: rationale and design of the randomized SORT OUT IX trial. *Am Heart J* 2019;213:1–7.
- [26] Okkels Jensen L. SORT OUT IX: A Randomized Trial Comparing a Polymer-free Coronary Drug-eluting Stent With an Ultra-thin Strut Bioreabsorbable Polymer-based Drug-eluting Stent in an All-comers Patient Population 2018. Available from <https://www.tctmd.com/slide/sort-out-ix-randomized-trial-comparing-polymer-free-coronary-drug-eluting-stent-ultra-thin>.
- [27] Lipinski MJ, Forrestal BJ, Iantorno M, Torguson R, Waksman RA. Comparison of the ultrathin Orsio hybrid sirolimus-eluting stent with contemporary drug-eluting stents: a meta-analysis of randomized controlled trials. *Cardiovasc Revasc Med* 2018;19(1 Pt A):5–11.
- [28] Iantorno M, Lipinski MJ, Garcia-Garcia HM, Forrestal BJ, Rogers T, Gajana D, et al. Meta-analysis of the impact of strut thickness on outcomes in patients with drug-eluting stents in a coronary artery. *Am J Cardiol* 2018;122(10):1652–60.
- [29] Bangalore S, Toklu B, Patel N, Feit F, Stone GW. Newer-generation ultrathin strut drug-eluting stents versus older second-generation thicker strut drug-eluting stents for coronary artery disease. *Circulation* 2018;138(20):2216–26.
- [30] Sergio B, Giovanna S, Bo L, Goran KO, Kristina H, Nils W, et al. Bioabsorbable polymer everolimus-eluting stents in patients with acute myocardial infarction: a report from the Swedish Coronary Angiography and Angioplasty Registry. *EuroIntervention* 2018;14(5):e562–e9.
- [31] Kobo O, Roguin A. Orsio: ultrathin bioabsorbable polymer sirolimus-eluting stent. *Future Cardiol* 2019;15(4):295–300.
- [32] Piccolo R, Heg D, Franzona A, Roffi M, Tuller D, Vuillienet A, et al. Biodegradable-polymer Sirolimus-eluting stents versus durable-polymer Everolimus-eluting stents in patients with acute ST-segment elevation myocardial infarction: insights from the 2-year follow-up of the BIOSCIENCE trial. *JACC Cardiovasc Interv* 2016;9(9):981–3.
- [33] Iglesias JF, Roffi M, Degrauwe S, Secco GG, Aminian A, Windecker S, et al. Orsio cobalt-chromium sirolimus-eluting stent: present and future perspectives. *Expert Rev Med Devices* 2017;14(10):773–88.
- [34] Roguin A, Kandzari DE, Marcusohn E, Koolen JJ, Doros G, Massaro JM, et al. Subgroup analysis comparing ultrathin, bioreabsorbable polymer sirolimus-eluting stents versus thin, durable polymer everolimus-eluting stents in acute coronary syndrome patients. *Circ Cardiovasc Interv* 2018;11(10):e007331.
- [35] Iglesias JF, Muller O, Zaugg S, Roffi M, Kurz DJ, Vuillienet A, et al. A comparison of an ultrathin-strut biodegradable polymer sirolimus-eluting stent with a durable polymer everolimus-eluting stent for patients with acute ST-segment elevation myocardial infarction undergoing primary percutaneous coronary intervention: rationale and design of the BIOSTEMI trial. *EuroIntervention* 2018;14(6):692–9.
- [36] Iglesias JF, Muller O, Heg D, Roffi M, Kurz DJ, Moarof I, et al. Biodegradable polymer sirolimus-eluting stents versus durable polymer everolimus-eluting stents in patients with ST-segment elevation myocardial infarction (BIOSTEMI): a single-blind, prospective, randomised superiority trial. *Lancet* 2019;394(10205):1243–53.
- [37] Moliterno DJ. The top papers of 2017: by subsequent citations and online views and downloads. *JACC Cardiovasc Interv* 2018;11(3):325–7.
- [38] Markovic S, Lutzner M, Dragomir S, Rottbauer W, Wöhrle J. Angiographic and clinical outcomes after recanalization of coronary chronic total occlusions with the Orsio sirolimus-eluting stent compared with the resolute zotarolimus-eluting stent. *Coron Artery Dis* 2017;28(5):376–80.
- [39] Sherbet DP, Christopoulos G, Karatasakis A, Daneke BA, Kotsia A, Navara R, et al. Optical coherence tomography findings after chronic total occlusion interventions: insights from the “Angiographic evaluation of the everolimus-eluting stent in chronic total occlusions” (ACE-CTO) study (NCT01012869). *Cardiovasc Revasc Med* 2016;17(7):444–9.
- [40] Teeuwen K, van der Schaaf RJ, Adriaenssens T, Koolen JJ, Smits PC, Henriques JP, et al. Randomized multicenter trial investigating angiographic outcomes of hybrid sirolimus-eluting stents with biodegradable polymer compared with everolimus-eluting stents with durable polymer in chronic total occlusions: the PRISON IV trial. *JACC Cardiovasc Interv* 2017;10(2):133–43.
- [41] Zivelonghi C, Agostoni P, Teeuwen K, van der Schaaf RJ, Henriques JPS, Vermeersch P, et al. 3-year clinical outcomes of the PRISON-IV trial: ultrathin struts versus conventional drug-eluting stents in total coronary occlusions. *JACC Cardiovasc Interv* 2019;12(17):1747–9.
- [42] Zivelonghi C, Teeuwen K, Agostoni P, van der Schaaf RJ, Ribichini F, Adriaenssens T, et al. Impact of ultra-thin struts on restenosis after chronic total occlusion recanalization: insights from the randomized PRISON IV trial. *J Interv Cardiol* 2018;31(5):580–7.
- [43] Teeuwen K, Spoormans EM, Bennett J, Dubois C, Desmet W, Ughi GJ, et al. Optical coherence tomography findings: insights from the “randomised multicentre trial investigating angiographic outcomes of hybrid sirolimus-eluting stents with biodegradable polymer compared with everolimus-eluting stents with durable polymer in chronic total occlusions” (PRISON IV) trial. *EuroIntervention* 2017;13(5):e522–e30.
- [44] Jensen LO, Maeng M, Raugaard B, Hansen KN, Kahlert J, Jensen SE, et al. Two-year outcome after biodegradable polymer sirolimus- and biolimus-eluting coronary stents (from the randomised SORT OUT VII trial). *EuroIntervention* 2018;13(13):1587–90.
- [45] Ellert J, Christiansen EH, Maeng M, Raugaard B, Jensen SE, Kristensen SD, et al. Impact of diabetes on clinical outcomes after revascularization with sirolimus-eluting and biolimus-eluting stents with biodegradable polymer from the SORT OUT VII trial. *Catheter Cardiovasc Interv* 2019;93(4):567–73.
- [46] Cassese S, Byrne RA, Tada T, Piniack S, Joner M, Ibrahim T, et al. Incidence and predictors of restenosis after coronary stenting in 10 004 patients with surveillance angiography. *Heart* 2014;100(2):153–9.
- [47] Iglesias JF, Heg D, Roffi M, Tuller D, Noble S, Muller O, et al. Long-term effect of ultrathin-strut versus thin-strut drug-eluting stents in patients with small vessel coronary artery disease undergoing percutaneous coronary intervention. *Circ Cardiovasc Interv* 2019;12(8):e008024.
- [48] Buiten RA, Ploumen EH, Zocca P, Doggen CJM, van der Heijden LC, Kok MM, et al. Outcomes in patients treated with thin-strut, very thin-strut, or ultrathin-strut drug-eluting stents in small coronary vessels: a prespecified analysis of the randomized BIO-RESORT trial. *JAMA Cardiol* 2019;4(7):659–69.
- [49] Jensen CJ, Richardt G, Tolg R, Erglis A, Skurk C, Jung W, et al. Angiographic and clinical performance of a paclitaxel-coated balloon compared to a second-generation sirolimus-eluting stent in patients with in-stent restenosis: the BIOLUX randomised controlled trial. *EuroIntervention* 2018;14(10):1096–103.
- [50] Neumann FJ, Sousa-Uva M, Ahlsson A, Alfonso F, Banning AP, Benedetto U, et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. *Kardiol Pol* 2018;76(12):1585–664.
- [51] Matkate S, Hanson M, Dallmann AC, Bentele M. Health economic evaluation of an ultrathin, bioreabsorbable-polymer sirolimus-eluting coronary stent compared to a thin, durable-polymer everolimus-eluting stent. *Cardiovasc Revasc Med* 2019;20(9):752–7.
- [52] Two-year results after implantation of an ultrathin, bioreabsorbable polymer sirolimus-eluting coronary stent compared with a thin, durable polymer everolimus-eluting stent: Health economic evaluation. In: Matkate S, Hanson M, Dallmann AC, Bentele M, editors. *American College of Cardiology (ACC 19). New Orleans, USA: Journal of American College of Cardiology (JACC); 2019 03/16/2019. [Available].*

- [53] Keehan SP, Stone DA, Poisal JA, Cuckler GA, Sisko AM, Smith SD, et al. National health expenditure projections, 2016–25: price increases, aging push sector to 20 percent of economy. *Health Aff (Millwood)* 2017;36(3):553–63.
- [54] Jakobsen L, Christiansen EH, Maeng M, Kristensen SD, Botker HE, Terkelsen CJ, et al. Randomized clinical comparison of the dual-therapy CD34 antibody-covered sirolimus-eluting combo stent with the sirolimus-eluting Orsiro stent in patients treated with percutaneous coronary intervention: rationale and study design of the Scandinavian Organization for Randomized Trials with clinical outcome (SORT OUT) X trial. *Am Heart J* 2018;202:49–53.